**REF** 70 380 Zn(x2)  
**REF** 70 402 TDA(x2)  
**REF** 70 542 JAMES(x2)

**REF** 70 422 VP 1(x2) + VP 2(x2)  
**REF** 70 442 NIT 1(x2) + NIT 2(x2)

**REF** 20 120 API 20 E™ REACTIFS  
API 20 E™ REAGENTS  
API 20 E™ REAGENZIEN  
API 20 E™ REACTIVOS  
API 20 E™ REATTIVI  
API 20 E™ REAGENTES  
API 20 E™ ΑΝΤΙΔΡΑΣΤΗΡΙΑ  
API 20 E™ REAGENSER  
API 20 E™ REAGENSERNE  
API 20 E™ ODCZYNNIKI

**REACTIFS**

oooooooo

**REAGENTS**

oooooooo

**REAGENZIEN**

oooooooo

**REACTIVOS**

oooooooo

**REATTIVI**

oooooooo

**REAGENTES**

oooooooo

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**COMPOSITION OF THE REAGENTS**

EN

<b>JAMES</b>	5 ml	R1: HCl 1N R2: Compound J 2183 (confidential)	100 ml 0.66 g
<b>NIT 1</b>	5 ml	Sulfanilic acid	0.4 g
<b>CORROSIVE:</b> <i>R34, S2, S23, S26</i>		Acetic acid H <sub>2</sub> O	30 g 70 ml
<b>NIT 2</b>	5 ml	N,N-dimethyl-1-naphthylamine	0.6 g
<b>CORROSIVE:</b> <i>R34, S2, S23, S26</i>		Acetic acid H <sub>2</sub> O	30 g 70 ml
<b>TDA</b>	5 ml	Ferric chloride (iron content) H <sub>2</sub> O	3.4 g 100 ml
<b>VP 1</b>	5 ml	Potassium hydroxide H <sub>2</sub> O	40 g 100 ml
<b>CORROSIVE:</b> <i>R35, S2, S26, S27, S37/39</i>			
<b>VP 2</b>	5 ml	$\alpha$ -naphthol Ethanol	6 g 100 ml
<b>IRRITANT / HIGHLY FLAMMABLE:</b> <i>R11, R36, S7, S16, S26</i>			
<b>Zn</b>	10 g	Zinc dust	
<b>HIGHLY FLAMMABLE / DANGEROUS FOR THE ENVIRONMENT:</b> <i>R15, R17, R50/53, S7/8, S43, S60, S61</i>			

The quantities indicated may be adjusted depending on the titer of the raw materials used.

*For more detailed information, consult the material safety data sheet available on request.*

**RISK AND SAFETY PHRASES**


EN

<b>R11</b>	Highly flammable.
<b>R15</b>	Contact with water liberates extremely flammable gases.
<b>R17</b>	Spontaneously flammable in air.
<b>R34</b>	Causes burns.
<b>R35</b>	Causes severe burns.
<b>R36</b>	Irritating to eyes.
<b>R50/53</b>	Very toxic to aquatic organisms, may cause long-term adverse effects in the aquatic environment.
<b>S2</b>	Keep out of the reach of children.
<b>S7</b>	Keep container tightly closed.
<b>S7/8</b>	Keep container tightly closed and dry.
<b>S16</b>	Keep away from sources of ignition - No smoking.
<b>S23</b>	Do not breathe vapors.
<b>S26</b>	In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.
<b>S27</b>	Take off immediately all contaminated clothing.
<b>S37/39</b>	Wear suitable gloves and eye/face protection.
<b>S43</b>	In case of fire involving metals, use powder extinguishers. Never use water.
<b>S60</b>	This product and its container must be disposed of as hazardous waste.
<b>S61</b>	Avoid release to the environment. Refer to safety data sheet.

## WARNINGS AND PRECAUTIONS

EN

- For *in vitro* diagnostic use and microbiological control.
- For professional use only.
- Do not allow reagents to come into contact with skin, eyes or clothing.
- Do not use reagents past the expiry date.
- Before use, check that the packaging and components are intact.
- Allow reagents to come to room temperature before use.
- Open ampules carefully as follows :

- Place the ampule in the ampule protector.
- Hold the protected ampule in one hand in a vertical position (white plastic cap uppermost).
- Press the cap down as far as possible.
- Position the thumb tip on the striated part of the cap and press forward to snap off the top of the ampule.
- Take the ampule out of the ampule protector and put the protector aside for subsequent use.

\* For ampule with no dropper-cap:

- Carefully remove the cap.

\* For ampule with dropper-cap:

- Turn the ampule upside down and maintain it in a vertical position.
- Squeeze gently on the cap to release a drop of reagent.

**NOTE:** For subsequent use, it is recommended to squeeze on the cap before turning the ampule upside down in order to suck back any spare drops of reagent and thus avoid sprinkling cap or fingers with the reagent.

## WARRANTY

bioMérieux disclaims all warranties, express or implied, including any implied warranties of MERCHANTABILITY AND FITNESS FOR A PARTICULAR USE. bioMérieux shall not be liable for any incidental or consequential damages. IN NO EVENT SHALL BIOMERIEUX'S LIABILITY TO CUSTOMER UNDER ANY CLAIM EXCEED A REFUND OF THE AMOUNT PAID TO BIOMERIEUX FOR THE PRODUCT OR SERVICE WHICH IS THE SUBJECT OF THE CLAIM.

## STORAGE CONDITIONS

EN

STORAGE CONDITIONS	In the dark	TDA / VP 1 / VP 2 / NIT 1 / NIT 2 / R1 and R2 of JAMES / API 20 E Reagent kit
STORAGE TEMPERATURE	2-8°C	VP 2 / NIT 2 / R2 of JAMES / API 20 E Reagent kit
	2-30°C	TDA / VP 1 / NIT 1 / R1 of JAMES
	8-30°C	Zn
STORAGE	Before opening	UNTIL THE EXPIRY DATE
	After opening (1) and transfer (except the non transferable reagents in the API 20 E Reagent kit)	<b>1 MONTH</b> (without exceeding the expiry date)
REAGENTS SENSITIVE TO LIGHT (2)		JAMES / VP 2

(1) **Record the date opened on the bottle or ampule label.**

(2) Check the appearance of the reagent before transferring it into the dropper-bottle. After transferring the contents of the ampule into the dropper-bottle, wrap the bottle in aluminium foil.

APPEARANCE OF REAGENTS	
JAMES	R1: Clear liquid – colorless R2: Before reconstitution: Whitish to yellowish After reconstitution: Yellow
NIT 1	Clear liquid – colorless
NIT 2	Clear liquid – colorless
TDA	Clear liquid – orangy yellow
VP 1	Clear liquid – colorless
VP 2	Clear liquid – light pink (*)

(\*) The shades of color may become deeper over time.

Important: Make sure that the reagents are put back in the refrigerator immediately after use.

## USE OF THE REAGENTS

Allow reagents to come to room temperature before use.

### 1. API 20 E™ reagent kit:

- Open the ampules of reagents as indicated in the paragraph "Warnings and Precautions" (ampule with dropper cap).

### 2. TDA, VP 1, VP 2, NIT 1, NIT 2:

- Open the ampule of reagent and transfer the contents into the dropper-bottle as indicated in the paragraph "Warnings and Precautions" (ampule with dropper cap).
- Dispense one drop of reagent.
- Carefully close the bottle after use and store it as indicated in the paragraph "Storage conditions".

**3. JAMES:**

- Open the ampule of solvent associated with the JAMES reagent (R1) as indicated in the paragraph "Warnings and Precautions" (ampule with no dropper-cap).
- Take up the contents of the ampule using a completely dry pipette and transfer this solvent into the dropper-bottle (R2).
- Fit the dropper to the bottle.
- Carefully close the bottle.
- Shake (as the bottle contains the dehydrated active ingredient).
- Wait approximately 10 minutes until the active ingredient is completely dissolved.
- Use the reagent thus reconstituted, carefully close the bottle and store it as indicated in the paragraph "Storage conditions".

**NOTE:** The JAMES reagent must only be used if it is pale yellow. If a pink color appears when the reagent is reconstituted with the solvent, wait until this pink color has completely disappeared before using the reagent.

**4. Zn:**

- Open the bottle.
- Take up an aliquot of powder (approximately 2-3 mg) using the spatula fixed to the cap and deposit this quantity in the reaction cupule.
- Carefully close the bottle after use and store it as indicated in the paragraph "Storage conditions".

**WASTE DISPOSAL**

Please refer to the paragraph "WASTE DISPOSAL" below the "INDEX OF SYMBOLS" table.

It is the responsibility of each laboratory to handle waste and effluents produced according to their type and degree of hazardousness and to treat and dispose of them (or have them treated and disposed of) in accordance with any applicable regulations.

***For more detailed information, consult the package insert of the strips with which the reagents are used.***

Printed in France

bioMérieux, the blue logo, API, API 20 E, RapiD 20 E and RAPIDEC are used, pending and/or registered trademarks belonging to bioMérieux SA or one of its subsidiaries.

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