TP N°1

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2.2 Cytological Technique

The characteristics of the transmission electron microscope impose severe restrictions on the nature and preparation of samples. Since electrons are easily absorbed by solid matter, only ultrathin sections (30 to 50 nm thick) can be examined. Additionally, since the major atoms of living matter (C, H, O, N) are transparent to electrons, it is also necessary to enhance the contrast of structures using heavy atoms that are opaque to electrons, which selectively bind to certain areas.

Steps of the protocol :

- **Fixation** : Traditional fixatives are not sufficient as they tend to destroy fine structures. Other chemicals, such as osmium tetroxide (OsO_4), glutaraldehyde, or potassium permanganate, must be used.
- **Dehydration** : The principle is the same as in classical histology, except that dehydration must end in a resin solvent used for embedding.
- **Embedding** : Embedding media are often resins of the araldite family. Solid and very hard blocks are obtained, including and impregnating the sample for analysis. The hardness of this material allows for extremely thin and regular sections.
- **Sectioning** : Ultrathin sections are produced using an ultramicrotome, which allows sections of 30 to 50 nm thick. These sections are collected on metal grids that allow electron passage.
- **Positive Contrast** : To enhance the contrast of electron-transparent samples, the ultrathin sections are impregnated with contrasting agents, usually heavy metal salts opaque to electrons, such as uranyl acetate or lead citrate, which preferentially bind to specific areas of the cell.

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Practical Session No. 1: Cell Study Techniques Using the Light Microscope

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- Identify the basic steps for preparing a temporary microscope slide.
- Use the light microscope properly to observe cells.
- Distinguish between plant and animal cells based on shape and components.

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Light microscope – glass slides – cover slips – forceps – dissecting needle – iodine solution – piece of onion – distilled water – absorbent paper or tissue.

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- Recall the concept of the cell as the structural and functional unit of living organisms.
- Show an image or model of a plant and an animal cell.
- Explain the importance of using the microscope to study cells.

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7. Take a small piece of onion.
8. Using forceps, separate a thin transparent layer from the inner epidermis.
9. Place it on a clean glass slide.
10. Add a drop of iodine solution (to stain the nucleus and cell wall).
11. Carefully place the cover slip to avoid air bubbles.
12. Wipe off the excess stain with absorbent paper.

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