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**Mohamed Khider University of Biskra**  
Faculty of Exact Sciences and Natural and Life Sciences  
Department of Natural and Life Sciences

Subject: General Microbiology  
Level: Second-year undergraduate  
Academic Year: 2024/2025  
Duration: 3 hours

**Lab Session No. 0: Precautions to Take During Microbiology Lab Practices**

### **I. Objectives**

- Understand the essential safety and hygiene rules in microbiology.
- Learn to prepare and maintain a sterile workstation.
- Familiarize students with the safe use of basic equipment (Bunsen burner, platinum loop, etc.).

### **II. General Recommendations**

*Personal attire and hygiene:*

- Wear a closed cotton lab coat.
- Tie back long hair and avoid loose clothing.
- Do not wear hand jewelry.
- Wash hands with soap and water or hand sanitizer before and after handling.

*Lab behavior:*

- No drinking, eating, or using phones.
- Strict adherence to schedules and teacher supervision.
- Immediately report any incident (broken glassware, spills, etc.).

*Workstation organization:*

- Disinfect the bench before and after each session with bleach or a bactericide.
- Dispose of biological waste in designated bins.
- Store equipment properly after use.

### **III. Equipment Overview and Usage**

*Workbench:*

- Clean with a sponge soaked in bleach.
  - Only keep necessary tools on the bench to avoid clutter.
- 

### **Bunsen Burner:**

#### *Before use:*

- Ensure the gas tap at the workstation (yellow) is off.
- Verify that the burner tap and air adjustment knob are closed.

#### *Safety measures:*

- Tie back long hair and avoid loose clothing to reduce fire risks.
- Keep flammable objects (alcohol bottle, crystallizer, etc.) away from the burner.

#### *Ignition:*

1. Open the gas supply:
  - Slowly open the gas tap at the workstation.
2. Light the flame:
  - Bring a lit match or lighter near the top of the burner.
  - Gradually open the burner tap until the flame ignites.
3. Adjust the flame:
  - Slowly turn the air intake knob to allow oxygen flow.
  - Create a blue flame (sterilizing flame) with a blue inner cone in the yellow flame center.

#### *Sterile zone use:*

- The sterile zone extends about 20 cm around the blue flame.
- Use this area to sterilize tools like the platinum loop by passing the tip through the flame until glowing.

#### *After use:*

- Turn off the burner tap first, followed by the gas tap.
- Ensure the burner cools before moving or storing it.

#### *Important precautions:*

- Never leave the burner unattended.
- Avoid working with flammable or loose clothing near the burner.
- In case of a problem (e.g., gas leak), turn off the burner immediately and alert the teacher.

-----

## **Platinum Loop:**

### *Description:*

- A platinum loop is essential in microbiology, consisting of three parts:
  - Handle: Held by hand, often made of metal or plastic.
  - Rod: Connects the handle to the loop.
  - Loop: Usually closed, 2 mm or larger in diameter, used for precise manipulations.
- Made of platinum, stainless steel, or plastic and can have various shapes (straight, L-shaped, or looped).

### *Main uses:*

- Transfer small amounts of liquid.
- Collect bacterial or fungal colonies.
- Spread a drop on a culture medium or microscope slide.

### *Procedure:*

#### **1. Sterilization before use:**

- Pass the loop through the flame of a Bunsen burner until it glows red-hot.
- Allow it to cool for a few seconds before contacting the sample to prevent heat damage.

#### **2. Handling:**

- Use the loop to collect or transfer the sample as needed.
- Always work within the sterile zone (about 20 cm around the blue flame).

#### **3. Sterilization after use:**

- Sterilize the loop again by passing it through the flame until glowing.
- Ensure all biological residues are destroyed to avoid cross-contamination.

-----

### *Important precautions:*

- Never place the hot loop on a non-sterile surface.
- Handle with care to avoid burns or contamination.
- Replace damaged or deformed loops immediately for accurate manipulations.